

WET SURFACES



# MANZIM SURGICO

Manufacturers & Exporters of Veterinary Instruments.



## ANIMAL HEALTH VETERINARY SURGICAL & DENTAL INSTRUMENTS CATALOGUE



# Company Profile

Manzim surgico was founded over a quarter century ago to fill the need for old world craftsmanship, combined with our German technology and production know-how to achieve this unique combination of skills and instruments of the finest quality Manzim Surgico has devised and created New process, New quality control System .This concept has been building an enviable reputation in operating rooms throughout the world.

## Products:

Manzim Surgico Strives to produce & Supply the latest and most comprehensive range of instruments, Due to rapid evolution in the field of biological science and constant revolution in surgical methodology,Manzim Surgico constantly review and update the product range with passage of time . All major production processes and performed in-house in production workshop.

## Activity:

Manzim Surgico manufacture and Supply top quality Re-useable type and Single-use Veterinary (Pet) Instruments.

## Quality Assurance:

Manzim Surgico prides itself on the rigorous quality control monitoring on every step from initial to final stage to ensure the supply of finest grade consistency instruments according to the required specifications supplied.

## Certified Quality Awards:

Manzim Surgico's Production and quality comply with MD EEC 92/93 ,ISO 9001,ISO 13485 ,FDA ,CGMP,CE MARK .

## Distribution:

Manzim Surgico's Instruments are distributed in South Korea,Spain, France, Germany, U.K, Rwanda, Australia,USA.

## Customer Services:

Manzim Surgico is whole-heartedly committed to total customer satisfaction as our policy is Maximum Customer Satisfaction through: Best Quality Instruments, Competitive prices and on time Deliveries.

**Contact No:** 0092-52-3253436

**Fax No:** 0092-52-3253436

**E-mail:** support@manzimsurgico.com

**info@manzimsurgico.com**

## **Order Procedure:**

When place an order by phone, fax or email please quote the following details.

- Order Number.
- Catalogue Number and description.
- Quantity
- Required Delivery Time.
- Special Marking (If any have)

## **Urgent Deliveries:**

In case of urgency, we can dispatch orders within 1 to 2 weeks time.

## **Terms And Conditions:**

The goods remain the property of Manzim Surgico until paid for in full. We only deal in terms of advance payment by T/T Telegraphic transfer or by irrevocable letter of credit L/C.

## **Guarantee:**

MANZIM SURGICO Guaranteed unconditionally for three years, they clear boil and autoclave testing and are also passivated. The stainless steel to use manufacture instruments is the most suitable for the intended use of them. The fact that it will resist corrosion as it is further tested by different kind of chemical and physical testing. We use 410ss series of steel for non-cutting products and 420ss and french ss for cutting and scaling products.

## **Facility Address:**

**M/S: MANZIM SURGICO.**

**KHATANA STREET NEAR ISLAM MASJID,  
P/O BOX 51310 S-I-E, SHATAB GARH, SIALKOT,  
PAKISTAN.**

**TEL:0092-52-3253436**

**FAX:0092-52-3253436**

**EMAIL:info@manzimsurgico.com**

**sales@manzimsurgi.com**

**Web:www.manzimsurgico.com**

**Contacting Person: M.MUHAMMAD IMRAN**

**(Chief Executive)**

**Mobile # 0092-300-5159397**

**0092-307-5311117**





**V-001**

**CANINE SPAY SET**

- A #3 Knife Handle
- B Mayo Scissors Straight 6-3/4"
- C Metzenbaum Scissors 7" Curved
- D Tissue Forceps 5-1/2", 1x2 Teeth
- E Brown Adson Forceps 4 3/4", 7x7 Teeth
- F Halsted Mosquito Forceps 5" Curved
- G Kelly Forceps 5-1/2" Curved
- H Roch-Carmalt Forceps 6-1/4" Cvd. + Str.
- I Backhaus Towel Clamp 3-1/2"
- J TC Olsen Hegar Needle Holder 6-1/2" TC
- K Allis Tissue Fcp 4x5 Teeth 6"
- L Snook Ovariectomy Hook

**V-002**

**FELINE SPAY SET**

- A #3 Knife Handle
- B Mayo Scissors Straight 6-3/4"
- C Metzenbaum Scissors 7" Curved
- D Tissue Forceps 5-1/2", 1x2 Teeth
- E Brown Adson Forceps 4 3/4", 7x7 Teeth
- F Halsted Mosquito Forceps 5" Curved
- G Kelly Forceps 5-1/2" Curved
- H Roch-Carmalt Forceps 6-1/4" Curved
- I Backhaus Towel Clamp 3-1/2"
- J TC Olsen Hegar Needle Holder 5-1/2" TC
- K Allis Tissue Fcp 4x5 Teeth 6"
- L Snook Ovariectomy Hook

**V-003**

**GENERAL SURGERY SET**

- A #3 Knife Handle
- B Mayo Scissors Straight 6-3/4"
- C Metzenbaum Scissors 7" Curved
- D Tissue Forceps 5-1/2", 1x2 Teeth
- E Brown Adson Forceps 4 3/4", 7x7 Teeth
- F Halsted Mosquito Forceps 5" Curved
- G Kelly Forceps 5-1/2" Curved
- H Roch-Carmalt Forceps 6-1/4" Curved
- I Backhaus Towel Clamp 3-1/2"
- J TC Olsen Hegar Needle Holder 5-1/2" TC
- K Snook Ovariectomy Hook
- L Allis Tissue Forceps 6", 5 X 6 Teeth
- M Sen Retractor
- N Suture Scissors, 8", sharp-Blunt





**V-004**

Exotic Animal Surgical Kit  
Kit includes:

- (1) Scissors, Surgical, 5.5"
- (1) Instrument Case, Sterilizing Boxes (with silicon Pads)50x152mm
- (1) Retractor, ALM Self-Retaining
- (1) Forceps, Dressing 12.5 cm straight
- (1) Scissors, Iris, 11 cm, straight tungsten carbide blades
- (1) Needle Holder, Olsen Hegar, with suture scissors
- (1) #3 Scalpel handle
- (1) Forceps, Micro Adson 12 cm
- (2) Forceps, Mosquito 3.5" curved
- (2) Forceps, Hartman Mosquito 3.5"



**V-005**

Mouse Vessel Cannulation Kit

- Instrument Case, Sterilizing Boxes (with silicon Pads)50x152mm
- Dumont Tweezers #5
- Vannas Scissors
- Alm Self-Retaining Retractor
- Fine, Serrated Forceps
- MicroFil, 28 gauge
- Vessel Cannulation Forceps, 0.5 mm Diameter





**V-006**

General Surgery Basic Kit  
Kit includes:

- (1) #3 Scalpel handle
- (1) Mayo scissors, 6.75", straight
- (1) Snook ovariectomy, 7"
- (3) Mosquito hemostatic forceps, 5", curved
- (2) Kelly hemostatic forceps, 5.5", straight
- (1) Metzenbaum scissors, curved, 7"
- (1) Brown-Adson forceps, 7x7 teeth
- (4) Rochester Carmalt forceps, 6.25", curved
- (4) Backhaus towel clamps, 5.25"
- (1) Olsen Hegar needle holders, 5.5"
- (1) Sterilization case



**MANZIM  
SURGICO**

**V-007**

Rat Dissecting Kit

- Dumont Tweezers #5
- Alm Self Retaining Retractor
- Vannas Scissors
- Needle Holder
- Iris Forceps, Curved, Serrated
- Blunt Probe, 1.0 mm Diameter
- SuperCut Iris Scissors, Straight, 10 cm





**V-008**

Feline Spay Pack  
Kit includes:

- (2) Kelly forceps, 5.5", curved
- (2) Kelly hemostatic forceps, 5.5", straight
- (2) Mosquito hemostatic forceps, 5", curved
- (2) Mosquito forceps, 5", straight
- (3) Jones towel clamp, 3.5"
- (1) Dressing forceps, 5.5"
- (1) #3 Scalpel handle
- (1) Adson Brown forceps, 4.75"
- (1) Mayo scissors, 6.75", straight
- (1) Mayo-Hegar needle holder, 6.7"
- (1) Olsen-Hegar needle holder, 6.5"
- (1) Snook ovariectomy, 7"



MANZIM  
SURGICO

**V-009**

Micro-surgical Kit  
Kit includes:

- (1) Tweezers, Dumont #5, 11cm non
- (1) Scissors, Spring, 14 cm curved
- (1) Needle Holder, Castroviejo 14 cm
- (1) Scissors, Tenotomy, 10 cm straight
- (1) Instrument Tray, autoclavable
- (2) Forceps, Micro Hemostatic, 12.5"
- (1) Knife Handle, Round Knurled
- (1) Blades, Miniature Edge 12-pack
- (1) Probe, 14 cm, stainless steel (no
- (1) Retractor, Weitlaner, 10.8 cm





**V-010**

Mouse Dissecting Kit

Dumont Tweezers #5

Wire Retractor

Vannas Scissors

Needle Holder

Iris Forceps, Curved, Serrated

Blunt Probe, 1.0 mm Diameter

Dissecting Scissors, 10 cm, Straight

Small Surgical Instrument Storage Portfolio



**MANZIM  
SURGICO**

**V-011**

Student Botanical Kit

1 Scissors, 14cm

1 Dropper

1 Scalpel Handle #3

4 Scalpel Blades

1 Teaser Needle, Angled

1 Teaser Needle, Straight

1 Mariam Tweezers, Angled

1 Tweezers, Straight

1 Instrument Storage Portfolio





**V-012**

Physiology II Kit

- Vannas Scissors
- Stevenson Retractor
- SuperCut Iris Scissors, Straight
- Iris Forceps, Curved, Serrated
- Rongeur 1.3 mm Jaw
- Adson Forceps, 1 x 2 Teeth
- Utility Scissors, Straight, Sharp/Blunt
- Olsen-Hegar Needle Holder
- Probe 1.0mm Diameter, Blunt
- Dumont Tweezers #5
- Operating Scissors, Straight, Sharp/Blunt
- Instrument Storage Portfolio





**V-013**

**Entomology Kit for Students**

**Kit Contains:**

- 2 Teaser Needles, Straight
- 2 Teaser Needles, Angled
- 1 Iris Scissors, 11.5 cm
- 1 Scalpel Handle #3
- 1 Adson Micro Forceps, 12 cm
- 1 Dropper
- 2 Micro Dissection Forceps, 12 cm, Straight
- 2 Micro Dissection Forceps, 12 cm, Curved
- 2 Entomology Forceps, 10 cm
- 1 Tweezers with Magnifying Mirror
- 1 Stainless Steel Ruler, 6"
- 10 Blades #10
- 1 Instrument Storage Portfolio



**MANZIM  
SURGICO**

**V-014**

**Dissection Kit for Students**

**Kit Contains:**

- 1 Scalpel Handle #3
- 1 Scalpel Handle #4
- 1 Dressing Forceps, 16 cm
- 1 Mayo Hager Needle Holder, 16 cm
- 1 Teaser Needle, Angled
- 1 Teaser Needle, Straight
- 1 Hooks with Chain
- 1 Tissue Forceps, 1:2, 16 cm
- 1 Dissecting Scissors, 16 cm, Sharp/Blunt
- 1 Kelly Forceps, 14 cm, Straight
- 1 Dropper
- 1 Stainless Steel Ruler, 6"
- 1 Grooved Director, Plain
- 1 Adson Tissue Forceps, 12 cm
- 4 Blades for Scalpel Handle #4
- 4 Blades for Scalpel Handle #3
- 1 Surgical Instrument Portfolio





**V-015**

Veterinary Periodontal Surgery Basic Kit  
Kit includes:

- (1) 503413 - Mirror, Cone socket
- (1) 503414 - Double end ex
- (1) 503415 - Periodontal pr
- (1) 500236 - #3 Scalpel har
- (1) 503417 - 1 Ochsenbein
- (1) 503419 - Kirkland gum l
- (1) 503420 - Orban scaler
- (1) 503421 - Jacquette scal
- (1) 503422 - Gracey curette
- (1) 503423 - Scaler, double
- (1) 503424 - Gracey curette
- (1) 503406 - Periosteal elev





**V-016**  
Student Veterinary Kit

Dissecting scissors, 4.5"  
Dissecting scissors, 4"  
Dressing forceps, 5.5"  
Micro dressing forceps, 4"  
#4 Knife handle  
Straight teasing needle  
Angled teasing needle  
#22 Surgical blade  
Canvas roll-up



**MANZIM  
SURGICO**



**V-017**  
Sterilization Mesh Cassette,  
4 x 4 x 2 cm (1.5 x 1.5 x .75")



**V-018**  
Sterilization Mesh Cassette,  
8 x 8 x 3.4 cm (3.1 x 3.1 x 1.4")




**V-019**  
Sterilization Mesh Cassette,  
10.5 x 7 x 2.5 cm (4.1 x 2.75 x 1")



V-020

CASTRATION EMASCULATORS INSTRUMENTS

- A WHITE Emasculator 12"
- B PLAIN Emasculator 12"
- C FRANK Emasculator 9"
- D Small WHITE Emasculator



V-021

STANDARD VETERINARY EYE SET

- A #3 Scalpel Handle
- B Wire Eye Speculum
- C Eye Dressing Forceps 4"
- D Chalzion Forceps
- E Stevens Tenotomy Scissors 4-1/8"
- F Castroviejo Needle Holder 5-1/2" Str.
- G Bishop-harman, Toothed Forceps
- H Baby Olsen-hegar N/Holder Tc Tip
- I Desmarres Retractor





V-022

STANDARD VETERINARY ORTHOPEDIC INSTRUMENT SET


- A Wire Cutter
- B Senn Retractor Blunt 6 3/8" (16.2cm), blunt prongs, double end
- C Senn Retractor Sharp 6 3/8" (16.2cm), sharp prongs, double end
- D Weitlaner Retractor 5 1/4" (13cm)
- E Freer Elevator 7" (17.8cm) double end, sharp and blunt blades
- F Bone Curette
- G Mallet
- H Bone File
- I Wire twister 6" (15.2) with Carb-N-Sert jaws
- J Kern Bone Holding Forceps Small 6 1/2" (15.2cm), with ratchet
- K Kern Bone Holding Forceps Large 9 1/2" (24.1cm), with ratchet
- L Bone Cutting Forceps
- M Pin Cutter
- N Osteotome





**V-023**
**SPECULA**


- A VIENNA Nasal Speculum, 5-3/4" (14.6cm), standard pattern, 1-1/8" (2.8cm) blades
- B VIENNA Nasal Speculum, 5-3/4" (14.6cm), standard pattern, 1-1/8" (2.8cm) blades (Vantage Brand)
- C VIENNA Nasal Speculum, 5-3/4" (14.6cm), light pattern, 1-1/8" (2.8cm) blades
- D KILLIAN Speculum, 6" (15.2cm), with set screw, 3-1/2" (8.9cm) blades
- E GRAVES Vaginal Speculum, 3-3/4" (14.6cm), small
- F GRAVES Vaginal Speculum, 3-3/4" (14.6cm), small (Vantage Brand)



MANZIM  
SURGICO

**V-024**
**ABDOMINAL PACK SET**

- A Balfour Retractor
- B Baby Balfour Retractor
- C Central Blade Only
- D De Bakey Forceps
- E Doyen Intestinal Clamps
- F Metzenbaum Scissor 18cm
- G T/c Olsen Hegar Needle Holder 18cm
- H Poole Suction Tube.






**V-025**

**NEUTER PACK SET**

- A Halstead Mosquito 12cm
- B Crile Wood N/holder Tc
- C Allis Tissue Forceps
- D Backhaus Towel Clips
- E Olsen-Hegar N/Holder
- F Metzenbaum Scissors
- G Adson Tissue Forceps
- H Scalpel Handle No.3 with scale



MANZIM  
SURGICO

**V-026**

**CRUCIATE PACK SET**

- A Baby Gelpi
- B Martin Cartilage Forceps
- C Deschamps Needle right small or large
- D Smilie Menisectomy Knife
- E Mini-Hohlmann
- F Hohmann Retractor





**V-027**

**CRUCIATE PACKS SET**

- |   |   |
|---|---|
| <p><b>A</b> Knife Handle, 10cm<br/> <b>B</b> Mayo Scissor 6.75" straight<br/> <b>C</b> Roger Wire Cutting Scissor 4.75"<br/> <b>D</b> Frazier Suction 10fr<br/> <b>E</b> Halsted-Mosquito Forceps 5" straight<br/> <b>F</b> Halsted Mosquito Forceps 5" curved<br/> <b>G</b> Crile Forceps 5.5" straight<br/> <b>H</b> Adson-Brown Forceps 4.75"<br/> <b>I</b> Olsen-Hegar Needle Holder 6.5"</p> | <p><b>J</b> MS-1193 Senn Retractor 6.25" (3 sharp prongs)<br/> <b>K</b> MS-1194 Mini Hohmann Retractor 6mm<br/> <b>L</b> MS-1195 Lempert Rongeur 6.5" curved<br/> <b>M</b> MS-1196 Wire Cutter T.C 14cm<br/> <b>N</b> MS-1197 Complete Hand Chuck Pin 3/32" smooth<br/> <b>O</b> MS-1198 Crimping Device<br/> <b>P</b> MS-1199 Tensioning Device<br/> <b>Q</b> MS-1200 Crimp Tubes<br/> <b>R</b> MS-1201 Cruciate Needles (set of 3 sizes) 80#<br/> <b>S</b> MS-1202 Monofilament Nylon Suture<br/> <b>T</b> MS-1203 Tray, Stainless Perforated</p> |
|---|---|





V-028

VETERINARY EXTRACTION BASIC SET-

- A #3 Scalpel Handle
- B Baby molt Periosteal Elevator
- C Extracting Forceps
- D Universal, mini, La Grange Scissors 4 1/4"
- E Olsen-Hegar Needle Holder Combined suture scissors 5 1/2"
- F Brown- Adson Tissue Forceps, 4 3/4"
- G Operating Scissors Straight 5 1/2"
- H Luxating Elevator Straight & Curved 2mm
- I Luxating Elevator Straight & Curved 3mm
- J Luxating Elevator Straight & Curved 5mm





V-029

**SPAYING & TEAT INSTRUMENTS**

- A COVAULT Ovariectomy Hook with probe tip, 8-1/4"
- B SNOOK Ovariectomy Hook, 8"
- C SNOOK Ovariectomy Hook, 8"
- D UDALL Teat Bistouri, 6-1/4", blunt
- E LICHTY Teat Knife, 5-1/2", sharp point
- F LICHTY Teat Knife, 5-1/2", blunt end
- G CORNELL Teat Curette, 6-1/4"
- H BOWER Obstetrical Forceps, 9-1/2", Straight, serrated jaws
- I BOWER Obstetrical Forceps, 9-1/2", curved, serrated jaws





**V-030**

**CANINE PERIODONTAL SET**

- A Mirror Handle Cone Socket
- B Mirror #4
- C Explorer 2
- D Probe 20-23 Explorer
- E Curette Columbia 13-14 Fine
- F Curette Columbia 4L-4R
- G Snook Ovariectomy Hook



**V-031**

**GROOMING INSTRUMENTS**

- A White Toenail Scissors 4-1/2"
- B Tissue & Cuticle Nipper 4" Convex Jaws
- C Tissue & Cuticle Nipper 4-1/2" Convex Jaws
- D Tissue & Cuticle Nipper 5" Convex Jaws
- E Tissue & Cuticle Nipper 5" Convex Hvy Pattern
- F Fetlock Clipper
- G Set of 2 Blades for Fetlock Clipper
- H Wing Nut for Fetlock Clipper






V-032

**CANINE & FELINE TARTER SCRAPERS**

- A Canine & Feline Tartar Scrapers Claw Type
- B Canine & Feline Tartar Scrapers Curved w Straight Edge
- C Canine & Feline Tartar Scrapers Angular, Pointed
- D Canine & Feline Tartar Scrapers Triangular
- E Canine & Feline Tartar Scrapers Curved w Angled Edge
- F Canine & Feline Tartar Scraper DE, Cvd w Angled Edges
- G Universal Scaler #12 - Double Ended, Solid
- H Towner Scaler #15 - Single Ended, Solid
- I Towner Scaler #15/33 - Double Ended, Solid
- J Universal Scaler #3 - Single Ended, Solid
- K Universal Scaler #33 - Single Ended, Solid
- L Universal Scaler #34 - Single Ended, Solid
- M Tartar Remover 4.5" Small
- N Tartar Remover 5.5" - Large






**V-033**

**RESTRAINT INSTRUMENTS**

- A Petite Mouth Gag 2-1/4" Lightweight
- B Feline Mouth Gag 3" Lightweight
- C Canine Mouth Gag 4" Small Size
- D Canine Mouth Gag 4" Medium Size
- E Canine Mouth Gag 6-3/4" Large Size



**V-034**

**VETERINARY PERIODONTAL BASIC KIT -MEDICAN**

- A Mirror Handle Complete
- B Probe/Explorer CP 12/23
- C Tower-Jaquette Scaler
- D Gracey Anterior Curette
- E Gracey Posterior Curette
- F Columbia Curette
- G Large Surgical Elevator
- H Small Surgical Elevator
- I Extra Surgical Elevator





CURETTES, SCALERS & PROBES



V-035 Jacquette Scaler #15/15 (33/33) Solid



V-036 Jacquette Scaler #25/35 (34/35) - Solid



V-037 Jacquette Scaler #2/3 (31/32) - Solid



V-038 Gracey Curette #11/12 - Solid



V-039 Gracey Curette #13/14 - Solid



V-040 Columbia Curette #13/14 - Solid Handle



V-041 McCall Curette #13/14 - Solid



V-042 Goldman-Fox Probe - Flat, Single Ended



V-043 Expro #23/Williams



V-044 Williams Probe



V-045 Williams/Fox Probe - Double Ended



V-046 Holmstrom Plugger/Spreader #90 - 45mm



V-047

**DENTAL PACK SET**

- A Winged Elevators
- B Dental Extracting Forceps
- C Root Elevator 76N
- D Explorer D/E
- E Curettes Columbia 13-14
- F Curettes Columbia 4L-4R
- G Scaler, Jacquette 2S-3S
- H Periosteal Elevator
- I T/c N/Holder





V-048

**RODENT PACK SET**


- A Mouth Gag
- B Check Retractors
- C Tooth Cutter
- D Tooth Extracting Forceps
- E Spatula
- F Molar Luxator 3mm
- G Incisor Luxator 5mm
- H Curettes





Veterinary Winged Dental Elevators

winged elevators are light-weight yet extremely durable. The sharp, tapered winged tip allows easy insertion into the PDL. The width of the winged tips have been shaped to match the size and shape of the tooth's root structure (sizes 1-6 mm). This unique winged shape increases control of the instrument and provides additional leverage to aid in extraction of the tooth.



V-049

Winged Elevators with Standard Handle

Winged Elevator 1mm



Winged Elevator 2mm



Winged Elevator 3mm



Winged Elevator 4mm



Winged Elevator 5mm



Winged Elevator 6mm



Straight Root Tip Pick

Used for assisting with extraction of retained root tips or root fragments.

Straight Root tip pick - stubby handle



Periosteal Elevator

Used for elevating the muco-periosteal flap during extractions on large breed dogs.

Periosteal Elevator - stubby handle





### Winged Elevators with Stubby Handle

The stubby handle elevators available as a 6 pc. set. See page 16

Winged Elevator 1mm



Winged Elevator 2mm



Winged Elevator 3mm



Winged Elevator 4mm



Winged Elevator 5mm




Winged Elevator 6mm



**V-050**

6 pc Stubby Handle Set includes stainless steel instrument holder/case





V-051 5pc Luxator instrument set (1-5mm Stubby)

1mm Luxator (Stubby Handle)



2mm Luxator (Stubby Handle)



3mm Luxator (Stubby Handle)



4mm Luxator (Stubby Handle)



5mm Luxator (Stubby Handle)



Luxators

Luxators, sizes 1-5mm are very fine and sharp. The 5 sizes allow their use for both feline and Canine teeth. They are used and recommended by dental specialist's for cutting the periodontal ligament and expanding the alveolus (tooth Socket) and are ideal for incisors and crowded teeth.

The stubby (shorter length and larger diameter) handle makes extraction so much safer as the users fingers can be placed at the instrument tip. Luxators are not elevators, we recommend winged elevators are used to elevate and extract the tooth.



**V-052**

**Explorer/Measuring Probe**

The probe end has periodontal markings etched in black. The fine but robust tip is ideal for periodontal probing in cats or dogs. The markings will not fade even after many autoclave cycles. The other end of the probe has a shepherd hook explorer to detect tooth defects both supra and subgingivally.



**V-053**

**Double Ended Periosteal Elevator**

The periosteal elevator is used to lift a muco-periosteal flap off the bone when performing surgical extractions. The fine double-ended 2/4 periosteal elevator is great for cats and small to medium dogs.



**V-054**

**Sickle Scaler**

This scaler is used to remove plaque and calculus from the tooth surface. The sickle scaler tip can be used to remove calculus from the developmental grooves and blood grooves in canine teeth.



**V-055**

**Jaquette Supragingival Scaler**

(Above the gum) The scaler is used to remove plaque and calculus from the tooth surface. It works well in tight occlusal spaces.



**V-056**

**Subgingival Curette**

(Below the gum) The curette is used to clean plaque and calculus from the tooth root surface.



**V-057**

**Minnesota Retractor**

Used to lift and protect the tissue and the user when performing surgical extractions whilst operating the high speed drill. One end holds the cheek back and the other holds the tongue down, a multi purpose retractor.







V-058

### Tartar Removing Forceps

The curved beak of the forceps allow for easy removal of heavy calculus from the tooth surface, decreasing the chance of fracturing the tooth.



V-059

### Extraction Forceps

Extraction forceps are used to grasp the tooth and remove from the socket after it is loosened. The small size and spring in the handle helps prevent excessive force and the risk of breaking the tooth. Overall length 125mm (5").



V-060

### Root Tip Extraction Forceps

Root tip extraction forceps have grooved tips with fine curved ends that meet when closed. They are ideal for root extraction or as an alternative to regular extraction forceps (D1012).  
6" / 150mm.



V-061

### Fragment Extraction Forceps

Root fragment forceps are ideal for extracting mobile broken root fragments and small, delicate teeth. The grooved and serrated jaws assist with grasping fragments.  
5 1/2" / 138mm.



Instruments Specially Designed for Felines and/or Very Small Dogs



**Straight Tip Elevators**

These super slim elevators are excellent for removing teeth in cats and small dogs. They are also useful in retrieving fractured roots. Tips can be sharpened using a dental stone.



**Elevator 1mm**

1mm wide straight tip elevator with stubby handle.

**Elevator 2mm**

2mm wide straight tip elevator with stubby handle.



**Probe/Explorer**

The probe end has periodontal markings etched in black. The fine but robust tip is ideal for periodontal probing in cats or dogs. The markings will not fade even after many autoclave cycles. The other end of the probe has a shepherd hook explorer to detect tooth defects both supra and subgingivally.



**Very Fine Periosteal Elevator**

Double ended elevator. The small delicate flat ends of this instrument make it ideal for raising gingival flaps in cats. The rounded tip helps prevent accidental slippage and tearing of the gingival flap.



**Resorptive Lesion Probe**

This probe can be used to detect resorptive lesions or residual calculus subgingivally. The probe can also be used to identify furcation involvement (loss of alveolar bone between tooth roots).



**Sickle Scaler**

This double ended hand instrument is used for the removal of supragingival (above the gum) plaque/calculus and stains.



**Curette**

This slimline curette is ideal for calculus and plaque removal both supra and subgingivally (above and below the gum). The fineness of the instrument is well suited for use in cats and small dogs.



**Extraction Forceps**

Extraction forceps - small. These fine extraction forceps are designed to grasp small teeth. The spring action allows their use without excessive force being applied to the tooth crown. 125 mm

**V-062**

Instruments mentioned on this page are available as a set .





### Curved Deciduous Elevators

The curved deciduous elevators are ideal for the removal of retained canines (milk teeth) and are available in three sizes suitable for the different sized deciduous teeth encountered.

The sharp elevator end cuts the PDL whilst the curve of the elevator tip follows the mesial curve and the other end the distal curve of the root. Curved deciduous elevators make the removal of retained teeth quick and easy because they follow the curve of the root.

**V-063**  
Curved Deciduous Elevator – Small (2mm)



**V-064**  
Curved Deciduous Elevator – Medium (3mm)



**V-065**  
Curved Deciduous Elevator – Large (4mm)



**V-066**  
Molt Periosteal Elevator – No. 9

The rounded edge of the periosteal instrument minimise the tearing of tissue when lifting the gum flap. The sharp end makes the initial lifting of the flap easier. Suitable for medium to large dogs.



**V-067**  
Small Dog or Cat Right Angle Extraction Forceps

Designed by IM3 for easier extraction of cat, small dog and deciduous teeth. Parallel right angle jaws reduce the risk of tooth fracture during extraction. By not meeting at the end, they reduce the risk of breaking off cats teeth. Small handles provide a more tactile feel.





### Instrument Sets

Instrument sets take the guess work out of choosing the correct dental instrument. Instruments are used throughout the world in teaching wetlabs and have been developed with the help of veterinary dental specialists. Made from high quality stainless steel, these instruments meet the needs of the modern day veterinarian.

#### V-068

##### 5 Piece Specialised Instrument Set

The set is intended as a supplement to ordinary tooth extraction sets. The curved deciduous elevators are ideal for the removal of retained canines (milk teeth). Each instrument is available individually.

The set consists of:

- Molt Periosteal Elevator - No 9
- Curved Deciduous Elevator - small
- Curved Deciduous Elevator - medium
- Curved Deciduous Elevator - large
- Extraction Forceps, right angled jaw



#### V-069

##### 14 Piece Instrument Set, Standard Handle Elevators

The set consists of:

- Winged elevators Size 1-5mm
  - Straight root tip pick
  - Tartar Removing Forceps
  - Extraction Forceps
  - Minnesota Retractor
  - Subgingival Curette
  - Jaquette Supraperiosteal Scaler
  - Sickle Scaler
  - Explorer / Measuring Probe
  - Double ended Periosteal Elevator
- Each instrument is available individually



#### V-070

##### 14 Piece Instrument Set, Stubby Handle Elevators

The set consists of:

- Winged elevators Size 1-5mm
  - Straight root tip pick
  - Tartar Removing Forceps
  - Extraction Forceps
  - Minnesota Retractor
  - Subgingival Curette
  - Jaquette Supraperiosteal Scaler
  - Sickle Scaler
  - Explorer / Measuring Probe
  - Double ended Periosteal Elevator
- Each instrument is available individually



#### V-071

##### Feline 8 Piece Instrument Set

The feline periodontal kit has been designed exclusively for use in cats and small dogs. Its design allows the practitioner to perform a complete oral examination (including detection of feline resorptive lesions) and cleaning of all tooth surfaces – both supra and subgingivally. All instruments are available individually.

The set consists of:

- 1 mm Wide Straight Tip Elevator
- 2 mm Wide Straight Tip Elevator
- Probe/Explorer
- Extraction Forceps
- Periosteal Elevator
- Resorptive Lesion Probe
- Sickle Scaler and Curette





**V-072**

**4 Piece Winged Elevator Set**

Winged Elevator sizes 1-4 mm in standard handle, supplied in a black pouch. Each elevator is available individually.

Standard Handle Set



**V-073**

**4 Piece Winged Elevator Set**

Winged Elevator sizes 1-4 mm in stubby handle, supplied in a black pouch. Each elevator is available individually.

Stubby Handle Set



**V-074**

**6 Piece Winged Elevator Set**

Winged Elevator sizes 1-6 mm in standard handle, supplied in a black pouch. Each elevator is available individually.

Standard Handle Set




**V-075**

**6 Piece Winged Elevator Set**

Winged Elevator sizes 1-6 mm stubby (shorter length and larger diameter) handle. The stubby handle makes extraction so much safer as the user's finger can be placed at the instrument tip. The new handle design is ideal for small or large hands. Supplied in a black pouch. Each elevator is available individually.

Stubby Handle Set






V-076

**Basic Dissecting Set**

A Collection of Instruments For Introductory Dissection vinyl pouch  
Scalpel 6" Scissors 4 1/2" Forceps student grade 4 1/2"  
Straight needle in Wooden Handle Bent needle in wooden handle  
Medicine Droppers.



V-077

**Dissecting Set**

Beginners For Elementary Junior High Or High School Life Science  
Vinyle pouch Scissors Forceps student /Scalpel with Replaceable  
blade Straight needle in wooden handle Bent needle in wooden  
handle Medicine dropper



V-078

**General Biology Dissecting Set**

Designed for work in general biology good Quality instruments  
scissor student Forceps fine point 4 1/2" S/S  
Scalpel 1 1/2 Blade ( with G pin S/S )  
Straight Needle with wooden handle Medicine Dropper  
single fold leatherette case.



V-079

**Fine Dissecting Set**

Designed for Freshman Courses in Biology  
Scissors Fine point 4 1/2" Forceps fine point with g pin 4 1/2"  
Scalpel replaceable Blade Teasing Needle Str  
Teasing Needle Cvd Medicine Dropper  
Single Fold Leatherette Case.



V-080

**Intermediate Dissecting Set**

Contains high quality instruments recommended for advance  
biology courses . Scissor surgical 4 1/2" Forceps fine point 4 1/2 "  
Surgeon scalpels (with # 22 blade) Teasing needle plastic handle Str  
Teasing needle Cvd Plastic handle Mall probe and seeker  
Medicine dropper Single fold leatherette case



V-081

**Tissue Culture Dissecting Set**

Features Stainless Steel Instrument in leatherette case  
Scissors M/Dissecting Iris 4 1/2 " Forceps watchmaker 4 1/4" Str  
Forceps watchmake 4 1/4" Cvd Scalpel micro dissecting 5"  
18x3 mm blade Teasing Needle Str S/S Teasing Needle Cvd S/S  
Single Fold leatherette case





**V-082**

**Anatomy Dissecting set**

A complete set Detailed Dissection for more advanced work in zoology and comparative anatomy .  
Surgical scissors 5 ½". Fine point scissors 4 ½ "  
Forceps fine point Cvd 4 ½ " Thumb forceps superior 5"  
Scalpel 1 ½" Blade .  
Cartilage knife/Flexible Probe. Triple chain Hook/Mall probe.  
Needle wooden handle Str. Needle wooden handle Str  
Medicine Dropper. Double fold leatherette case



**V-083**

**Graduate Dissecting Set**

Design for the biology graduate students  
Surgical scissors 5 ½". Fine point scissors 4 ½".  
Forceps fine point Str 4 ½" Forceps fine point Cvd 4 ½ "  
Tripple chain hook. cartilage knife scalpel stain finish 1 ½" blade  
forceps cartilage 5" 2 needle holder with straight needle .  
Flexible probes. Probes with hook. Probe & seeker  
Dropping Pipe Double Fold Leatherette case



**V-084**

**Micro Dissecting Set**

Feature stainless steel Instruments in Leatherette case.  
Scissors M/Dissecting Iris 4 ½" Scissors sharp point Iris 4 ½ "  
Forceps Micro Dissecting 4 1/4" Straight Forceps micro dissecting 4"  
Curved Scalpel micro dissecting 18 x 3 mm blade  
scalpel micro dissecting 30 x 3 mm blade Dissecting needle Str .  
Dissecting needle Cvd. Dissecting needle Full-Spear point  
Dissecting needle Half-Spear point In double folding case.




**V-085**

**Large Animal set**

This set provide all the Necessary instruments Dissection of large Animal such as Cat, Shark and Rabbits Scissor surgical 5 ½ "  
D/Scissors superior 4 ½ " Scissors fine point Str 4 ½ "  
Forceps large Tissue 1x2 teeth 6" Forceps utility Cvd point 6" .  
Forceps fine point str 4 ½" Scalpel screw Lock with # 22 blade .  
Knife cartilage 6 ½ ". Probe sharp point angular. Section lifter S/S  
Teasing needle Str Wooden handle. Teasing Needle Cvd Wooden handle  
Medicine dropper.



**V-086**  
**RANGE OF SPATULAS**  
 7 parts filled, in a case with press stud  
 Content  
 1 straight blunt tweezers, 145 mm  
 1 micro double spatula, 150/4 mm  
 1 double spatula, 150/9 mm  
 1 micro-powder spatula, 150/4 mm  
 1 powder spatula, 150/9 mm  
 1 micro-spoon spatula, 150/5 mm  
 1 micro-spoon spatula, 150/7 mm



**V-087**  
**Jewellery Forceps Set**



**V-088**  
**DISSECTING SET**  
 6 parts or 12 parts filled, in pouch  
 Content  
 1 powder spatula, 150/9 mm  
 1 straight scissors, 100 mm  
 1 spoon spatula, 150 mm  
 1 scalpel 125 mm  
 1 flexible section lifter, 150 mm  
 1 pointed tweezers, 130 mm



**V-089**  
**Needle Set**

**V-090**  
**BOWMAN LACRIMAL PROBES**  
 Double Ended, Sterling Silver



Bowman Probe, 000-0000  
 Bowman Probe, 00-0  
 Bowman Probe, 1-2  
 Bowman Probe, 3-4  
 Bowman Probe, 5-6  
 Bowman Probe, 7-8

**V-091**  
**WILLIAM LACRIMAL PROBES**  
 Double Ended, Sterling Silver



William Probe, 000-0000  
 William Probe, 00-0  
 William Probe, 1-2  
 William Probe, 3-4  
 William Probe, 5-6  
 William Probe, 7-8



V-092  
Fish Aquarium Set




V-093  
Fish Aquarium Set



V-094  
Fish Aquarium Set



V-095  
Fish Aquarium Set



V-096  
Fish Aquarium Set



V-097  
Fish Aquarium Set





V-098



V-099



V-100



V-101



V-102  
Dissection Premium Kit



**V-103**

**DRESSING SET - 8 PIECES**

- Kidney Dish 20 cm (8") - medium
- Gallipot 150 ml - 6 oz - medium
- Kocher artery forceps - straight - 16 cm
- Dressing Forceps (Thumb) - 16 cm
- English TOE dissecting forceps - plain - 14 cm
- Mayo scissors - straight - 15 cm
- Rampléy sponge holding forceps - str. - 18 cm
- Instrument tray with lid 25x30x6.2 cm



**V-104**

**SUTURE SET - 12 PIECES**

- Kidney Dish 20 cm (8") - medium
- Gallipot 150 ml - 6 oz - medium
- Mayo Hegar needle holder - 18 cm
- Kocher artery forceps - straight - 16 cm
- Kocher artery forceps curved - 16 cm
- Dressing Forceps (Thumb) - 16 cm
- English TOE dissecting forceps - plain - 14 cm
- Dissecting forceps toothed 1x2 - 14 cm
- Suture scissors - Heaths Ligature - 15 cm
- Mayo scissors - straight - 15 cm
- Rampléy sponge holding forceps - str. - 18 cm
- Instrument tray with lid 30x20x5 cm



**V-105**

**MINOR DELIVERY SET - 25 PIECES**

- Kidney Dish 20 cm (8") - medium
- Kidney Dish 25 cm (10") - large
- Gallipot 150 ml - 6 oz - medium
- Brown Sadler episiotomy scissors - 14 cm
- Kocher artery forceps - straight - 16 cm
- Kocher artery forceps straight - 18 cm
- Kocher artery forceps curved - 16 cm
- English TOE dissecting forceps - plain - 14 cm
- English TOE dissecting forceps - plain - 18 cm
- Dissecting forceps toothed 1x2 - 14 cm
- Operating scissors sharp/blunt - straight - 15 cm 2 pcs.
- Operating scissors sharp/blunt - curved - 15 cm 2 pcs.
- Rampléy sponge holding forceps - str - 24 cm 2 pcs.
- Museux Vulsellum forceps 8 mm - str. - 24 cm
- Cusco vaginal specula - small
- Cusco vaginal specula - medium
- Cusco vaginal specula - large
- Sims vaginal specula - small
- Sims vaginal specula - medium
- Sims vaginal specula - large
- Umbilical scissors USA Patent - 10.5 cm
- Instrument tray with lid 45x30x6.2 cm





**V-106**  
 PROBE  
 BUTTERFLY PROBE  
 STRAIGHT BLADE  
 POT-BELLIED BLADE  
 ANATOMY FORCEPS  
 SURGERY FORCEPS  
 HOOK SHARP  
 FORCEPS THIN END  
 SCISSORS BLUNT/SHARP



**V-107**  
 PROBE  
 BUTTERFLY PROBE  
 STRAIGHT BLADE  
 POT-BELLIED BLADE  
 ANATOMY FORCEPS  
 SURGERY FORCEPS  
 KLEMMER FORCEPS  
 STITCH REMOVER  
 FORCEPS THIN END  
 SCISSORS BLUNT/SHARP  
 MATHIEU NEEDLE HOLDER

MANZIM  
 SURGICO



**V-108**  
 PROBE  
 BUTTERFLY PROBE  
 STRAIGHT BLADE  
 POT-BELLIED BLADE  
 ANATOMY FORCEPS  
 SURGERY FORCEPS  
 STITCH REMOVER  
 FORCEPS THIN END  
 SCISSORS BLUNT/SHARP  
 MAYO-HEGAR HOLDER



**V-109**  
 PROBE  
 BUTTERFLY PROBE  
 STRAIGHT BLADE  
 POT-BELLIED BLADE  
 ANATOMY FORCEPS  
 SURGERY FORCEPS  
 KLEMMER FORCEPS  
 STITCH REMOVER  
 FORCEPS THIN END  
 SCISSORS BLUNT/SHARP  
 MATHIEU NEEDLE HOLDER





**V-110 Set of chrome plated Hegar Dilators**  
Include 14 dilators Ø 4-17 mm

Sims hystrometer - malleable -  
33 cm - silver

Bosch/Sims scissors - 23 cm

Schroeder forceps - 25 cm  
(Pozzi-Palmer)

Foerster forceps - 25 cm

**GYNAECOLOGY KIT**  
Includes 4 instruments marked  
with in a nylon bag.





# ***MANZIM SURGICO***

**We Are Manufacturer & Exporter of Veterinary (Pet) Instruments.**

**ADDRESS :** R/CN 394 Near Islam Masjid Shatab Garh  
P/o Small Industrial Estate Sialkot. 51310 Pakistan

**TEL :** 0092-52-3253436

**FAX :** 0092-52-3253436

**MOBILE#** 0092-300-5159397  
0092-307-5311117

**EMAIL :** [info@manzimsurgico.com](mailto:info@manzimsurgico.com)  
[sales@manzimsurgico.com](mailto:sales@manzimsurgico.com)

**WEB:** [www.manzimsurgico.com](http://www.manzimsurgico.com)