



VETERINARY INSTRUMENTS CATALOGUE



www.manzimsurgico.com

we are manufacturer & exporter of veterinary instruments

E-mail: info@manzimsurgico.com

Company Profile

Manzim surgico was founded over a quarter century ago to fill the need for old world craftsmanship, combined with our German technology and production know-how to achieve this unique combination of skills and instruments of the finest quality ManzimSurgico has devised and created New process, New quality control System .This concept has been building an enviable reputation in operating rooms throughout the world.

Products:

Manzim Surgico Strives to produce & Supply the latest and most comprehensive range of instruments, Due to rapid evolution in the field of biological science and constant revolution in surgical methodology,Manzim Surgico constantly review and update the product range with passage of time. All major production processes and performed in-house in production workshop.

Activity:

Manzim Surgico manufacture and Supply top quality Re-useable type Veterinary instruments .

Quality Assurance:

Manzim Surgico prides itself on the rigorous quality control monitoring on every step from initial to final stage to ensure the supply of finest grade consistency instruments according to the required specifications supplied.

Certified Quality Awards:

Manzim Surgico's Production and quality comply with MD EEC 92/93 ,ISO 9001,ISO 13485 ,FDA ,CGMP,CE MARK .

Distribution:

Manzim Surgico's Instruments are distributed in South Korea,Spain, France, Germany, U.K, Rwanda, Australia,USA.

Customer Services:

Manzim Surgico is whole-heartedly committed to total customer satisfaction as our policy is Maximum Customer Satisfaction through: Best Quality Instruments, Competitive prices and on time Deliveries.

Contact No: 0092-52-3253436

Fax No: 0092-52-3253436

E-mail: support@manzimsurgico.com

info@manzimsurgico.com

Order Procedure:

When place an order by phone, fax or email please quote the following details.

- Order Number.
- Catalogue Number and description.
- Quantity
- Required Delivery Time.
- Special Marking (If any have)

Urgent Delivieries:

In case of urgency, we can dispatch orders within 1 to 2 weeks time.

Terms And Conditions:

The goods remain the property of Manzim Surgico until paid for in full. We only deal in terms of advance payment by T/T Telegraphic transfer or by irrevocable letter of credit L/C.

Gurantee:

MANZIM SURGICO Guaranteed unconditionally for three years, they clear boil and autoclave testing and are also passivated .The stainless steel to use manufacture instruments is the most suitable for the intended use of them. The fact that it will resist corrosion as it is further tested by different kind of chemical and physical testing. We use 410ss series of steel for non-cutting products and 420ss for cutting products.

Facility Address:

M/S: MANZIM SURGICO.

**KHATANA STREET NEAR ISLAM MASJID,
P/O BOX 51310 S-I-E, SHATAB GARH, SIALKOT,
PAKISTAN.**

TEL:0092-52-3253436

FAX:0092-52-3253436

EMAIL:info@manzimsurgico.com

sales@manzimsurgi.com

Web:www.manzimsurgico.com

**Contacting Person: M.MUHAMMAD IMRAN
(Chief Executive)**

**Mobile # 0092-300-5159397
0092-307-5311117**


INST CASTRATION ESAMULATORS



Bloodless Castrator For Lambs

MS-001-01

**Bloodless Castrator for Lambs
9"**



Bloodless Castrators For Bulls

MS-001-02

**Castration
18"**



Castrator

MS-001-03

**Castration 9",
12", 16", 19"**



Castrators

MS-001 - 04

**Castration
9"**



Cox Bloodless Castrator

MS-001-05

**Castration
16"**



New Barry Castrating Knife

MS-001-06

**Castration
10.5"**

INST CASTRATION EMASCULATORS




Inst Castration Emasculators
MS-001-07
Inst Castration
Emasculators With wooden
handle
30 cm, 113/4



Inst Castration Emasculator
MS-001-08
Inst Castration
Emasculator 23 cm, 9"



Castrating Forceps
MS-001-09
Castrating
forceps 35 cm,
14"



Verboczy
MS-001-10
With clamping parts 2
27 cm, 102/8



Haussmann
MS-001-11
Straight
27 cm, 105/8, 31 cm,
121/4



Haussmann
MS-001-12
Curved
19 cm, 71/2, 27 cm,
105/8

KNIFES



**F Dick Ascot Hoof Knife
MS-001-13
Knives**



**Left-Handed Hoof Knife
With Wooden Handle
MS-001-14
Knives**



**Right-Handed Hoof Knife
With Wooden Handle
MS-001-15
Right-Handed Hoof Knife With
Wooden Handle**



**Knife Single Edge Plastic Handle
MS-001-16
Knife Single Edge Plastic
Handle**



**Hoof Knife W Stainless Blade - Right
Hand
MS-001-17
Hoof Knife With Stainless Blade -
Right Hand**



**Cottam Farriers Hoof Knife
MS-001-18
Cottam Farriers Hoof Knife**

KNIFES



**Oval Hoof Knife
MS-001-19**



**Hoof Pick
MS-001-20**



**Hoof Knife
MS-001-21**



**With Wooden Handle Left & Right
Side**

MS-001-22

**With wooden handle left & right
side 19 cm , 7 1/2**



With Wooden Handle

MS-001-23

**With wooden
handle 20 cm, 8"**



**Hoof Knife
MS-001-24
Stainless steel
handle 20 cm, 8"**

HOOF NIPPERS



Half Round Hoof Nipper
MS-001-25

The smallest of the nippers. Ideal for foals and any of the smaller-hoofed horses
12"



Hoof Nipper
MS-001-26

The 14" nipper is considered the workhorse of the nippers. It is especially useful when working on hard hooves.
Designed to be used on all breeds of horses
14"



Hoof Nipper
MS-001-27

The 15" nipper is the largest of nippers. Farriers doing larger breeds and large amounts of trimming consider this nipper to be their most important tool
15cm



Easy Racetrack Hoof Nipper
MS-001-28

The 14" race track easy Nipper combines the comfort of the easy handle with better balance.
14"



Easy Hoof Nipper
MS-001-29

The 15" easy Nipper is the largest of the 'Easy' range that combines the comfort of the easy handle with better balance, increased leverage.
15"



Semi-Polished Hoof Nipper
MS-001-30

The 12" is the smallest of the nippers and is ideal for smaller-hoofed horses.
12"

HOOF NIPPERS



Racetrack Semi-Polished Hoof Nipper

MS-001-31

The 14" race track nipper forged from high quality tool steel, machined to precise tolerances and heat treated for wear resistance.

14"



Polished Hoof Nipper

MS-001-32

The 12" is the smallest of the nippers and is ideal for smaller-hoofed horses.

12"



Racetrack Polished Hoof Nipper

MS-001-33

The 14" race track nipper forged from high quality tool steel, machined to precise tolerances and heat treated for wear resistance.

14"



Plastic Grip Hoof Nipper

MS-001-34

Designed exclusively for trimming hooves. This 14-inch long professional nipper features a wide 1" flat head and thin heat-treated jaw. Equipped with cushion-grip handles.

14"



Polished Hoof Nipper

MS-001-35

The Nordic Forge 15" Nipper is forged from high quality tool steel, machined to precise tolerances and heat treated for wear resistance. Offers a very sharp cut.

15"



Hoof Nipper

MS-001-36

14"

HOOF NIPPERS



Nail Tong
MS-001-37
14"



Apex 14DN 14-Inch Hoof Nippers
MS-001-38
Apex 14DN 14-inch Hoof
Nippers14"



Diamond Cooper Hand Tools FN15
15-In Farrier Nipper
MS-001-39



Hoof Nipper
MS-001-40



Hoof Nipper
MS-001-41
12"



One Handed 12 Nipper
MS-001-42
One Handed 12 Nipper

HOOF NIPPERS



Riveting Tongs
MS-001-43
Riveting Tongs
14"



Alligator Clinchers
MS-001-44
Alligator
Clinchers12"



Alligator Tongs
MS-001-45
Alligator
Tongs14"

**MANZIM
SURGICO**

HOOF TESTERS



Hoof Tester
MS-001-46
12"



Hoof Tester
MS-001-47
13"



Hoof Tester
MS-001-48
14"



MANZIM
SURGICO

FARRIERS HAMMERS



**Apex FH14 14 Oz Farrier Driving
Hammer
MS-001-49
Farrier Hammer**



**Diamond Farrier Driving Ham
MS-001-50
Farrier Hammer**



**Driving Hammer
MS-001-51
Farriers Hammer
Driving Hammer
14"**



**Farrier Hammer
MS-001-52
Farrier Hammer
12"**

**MANZIM
SURGICO**

BARNES DEHORNER



**Barnes Dehorner
MS-001-53
Barnes Dehorner
17"**



**Barnes Dehorner
MS-001-54
Barnes Dehorner
20"**



**Barnes Dehorner
MS-001-55
Made of solid copper with 20 steel
handle. Set includes 5/8", 3/4", and
7/8".**



**Barnes Dehorner Wooden Handle
MS-001-56
Barnes Dehorner
13"**



**Convex Dehorner
MS-001-57
Barnes Dehorner
25"**



**Convex Junior Dehorner
MS-001-58
Barnes Dehorner
14"**

BARNES DEHORNER



**Dehorner
MS-001-59
Barnes Dehorner
17"**



**Improved Keystone Dehorner
MS-001-60
Barnes
Dehorner24"**



**Keystone Dehorner
MS-001-61
Barnes
Dehorner44.5"**



**Knuckle Saver Precision Barnes
Dehorner
MS-001-62
Barnes Dehorner
17"**



**Knuckle Saver Precision Barnes
Dehorner
MS-001-63
Barnes
Dehorner20"**



**Professional Barnes Style Dehorner
MS-001-64
Barnes
Dehorner8",
11", 14"**

BULL HOLDERS



**Brass Nose Leader
MS-001-65
Bull Holder**



**Nose Leader
MS-001-66
Bull Holder**



**Bull Leader
MS-001-67
Bull Holder**



**Spring Nose Leader
MS-001-68
Bull Holder**



**Heavy Duty Bull Leader
MS-001-69
Bull Holder**



**Nose Leader
MS-001-70
Bull
Holder 7"**

BULL HOLDERS



Bull Holder With Chain
MS-001-71



Bull Ring Holder
MS-001-72



Bull Nose Tweezers
MS-001-73



**MANZIM
SURGICO**

BULL RINGS



Bull Ring - Brass - Small 2 1/2" X 5/16"

MS-001-74

Bull Ring - Brass - Small 2 1/2" X 5/16"



Bull Ring - Brass - Medium 3" X 3/8"

MS-001-75

Bull Ring - Brass - Medium 3" X 3/8"



Bull Ring - Brass - Dairy 3" X 5/16"

MS-001-76

Bull Ring - Brass - Dairy 3" X 5/16"



Bull Ring - Stainless - 2 1/2"

MS-001-77

Bull Ring - Stainless - 2 1/2"



Bull Ring - Stainless - 2 3/4"

MS-001-78

Bull Ring - Stainless - 2 3/4"



Bull Ring - Stainless - 3"


MS-001-79

Bull Ring - Stainless - 3"

BULL RINGS



Copper Bull Ring 2.25" , 2.75" , 3"
MS-001-80



**MANZIM
SURGICO**

EAR NOTCHER



Ear Notcher, Small 'V', (5/16" At Base)

MS-001-81
Ear Notcher



Ear Notcher V Large

MS-001-82
Ear notcher v shape
6", 7", 8", 9"



Ear Notcher U Shape

MS-001-83
Ear notcher
U shape stainless steel.
6", 7", 8", 9"



Ear Notcher - U Notch Large

MS-001-84

Ear Notcher - U Notch Large



Medium & Large Ear Notcher

MS-001-85

Medium & Large Ear Notcher
with Inverted U-Shaped Die

EMESCULATOR



**Sand S Castration Clamp
MS-001-86
Emesculator
12.5"**



**Reimer S Emasculator
MS-001-87
Emesculator
13"**



**Serra Emasculator
MS-001-88
Emesculator
14"**



**Verboeczy Emasculator
MS-001-89
11"**



**Veterinary Emasculator
MS-001-90
Emesculator
12"**

RING APPLICATORS



Bull Nose Punch
MS-001-91
Bull Nose Punch



Bull Ring Applicator
MS-001-92
10"



Castration Ring Applicator
MS-001-93



Elastrator Standard
MS-001-94




Model Bison, Ring Applicator
MS-001-95



Pig Nose Ring Applicators
MS-001-96

RING APPLICATORS



Castration Ring Applicator

MS-001-97

**Castration Ring
Applicator 9"**

**MANZIM
SURGICO**

MILK SUCKING PREVENTER




**Milk Sucking Preventer
MS-001-98**



**Milk Sucking Preventer
MS-001-99**



**Plastic Anti Suckler
MS-001-100**



**MANZIM
SURGICO**

EQUESTRIA



Equestrian Equipment
MS-002-01

Equestrian Equipment
MS-002-02
Equestrian
Equipment 22 cm, 8
3/4

Equestrian Equipment
MS-002-03
22 cm, 8 3/4




Equestrian Equipment
MS-002-03
22 cm, 8 3/4

Equestrian Equipment
MS-002-04
22 cm, 8 3/4

Equestrian Equipment
MS-002-05
22 cm, 8 3/4

EQUESTRIA



Equestrian Equipment
MS-002-06

22 cm, 8 3/4



Equestrian Equipment
MS-002-07

22 cm, 8 3/4

**MANZIM
SURGICO**

DEHORNING INSTRUMENTS



**Caw Cutter, Side Action
MS-002-08 Dehorning
Instrument
60 cm, 231/2**



**Barnes Dehorner (American Pattern)
MS-002-09
Dehorning Instrument
60 cm, 231/2**



**For Claves Age 3 - 6 Months
MS-002-10**



**Hoof- And Claw Cutter
MS-002-11
Dehorning Instrument
41cm,16**



**Krause
MS-002-12
Milk Sucking
Preventer 41 cm, 16**



**Wide Toothed
MS-002-13
Clipper with wooden
handle 41 cm, 16**

BALLING GUNS



Standard Balling Gun
MS-002-14

(Cotran) 15.5" long, holds up to 1 ounce size capsule. Metal handle with spring.



Balling Gun - Calf Size
MS-002-15

(Stone Manufacturing Company) New non-injurious bolus dispenser with flexible plastic head.



Balling Gun - Medium
MS-002-16

New non-injurious bolus dispenser with flexible plastic head.



Farmers Balling Gun
MS-002-17

Single dose balling gun with 1" head features heavily chrome-plated metal to resist corrosion for long, trouble-free service and reinforced heads for extra strength. 18" balling gun includes bolus spring clips.



18"X1/2" Balling Gun W/ Chrome Plated Head
MS-002-18

Single dose balling gun with 1/2" head features heavily chrome-plated metal to resist corrosion for long, trouble-free service and reinforced heads for extra strength. 18" balling gun includes bolus spring clips.



18"X5/8" Balling Gun W/ Chrome Plated Head
MS-002-19

Single dose balling gun with 5/8" head features heavily chrome-plated metal to resist corrosion for long, trouble-free service and reinforced heads for extra strength. 18" balling gun includes bolus spring clips.

BALLING GUNS



**18"X7/16" Balling Gun W/ Chrome Plated Head
MS-002-20**

Single dose balling gun with 7/16" head features heavily chrome-plated metal to resist corrosion for long, trouble-free service and reinforced heads for extra strength. 18" balling gun includes bolus spring clips.



**18" Balling Gun Frame Only
MS-002-21**

The WIMCO Cler-Head Balling Gun frame accepts interchangeable plastic heads up to 1" in 1/16" increments to fit boluses up to 480 grams and many common rumen magnets. Order Frames and Heads separately.



**Dose Syringe With Straight Pipe
MS-002-22**

6 oz. dose syringe with 9" straight pipe. Used to dispense liquid medication orally. Plunger rod is calibrated in 1/2 oz. increments. Heavy-duty, nickel finish brass construction. Easy disassembly for cleaning. Leak proof operation. Metal plunger with rubber O-ring.

COW PULLER



Ratch-A-Pull Dual Action
MS-002-23

New, improved model double action puller. Refined milled grooves on shaft with no sharp edges. Wide cast aluminum breech spanner and longer solid shaft for those big calves. Puller disassembles for storage. Heavy-duty construction throughout. Order OB chains separately.



OB Handle
MS-002-24

Hook, nickel finish, is designed to grasp chain firmly at any length.



OB T-Bar Handle
MS-002-25

Easy grip two-handed pull design provides extra leverage. Hook will not spring or let chain slip.



Wire Saw
MS-002-26

Abrasive wire saw for rapid cutting. Wire saw can be useful when dehorning or for removal of dead fetus from the womb. Handles sold separately.

SHEEP SHEARS



Burgon & Ball Double Bow Bent Shear

MS-002-27
Sheep Shear
5.5"



Double Bow & Single Bow Sheep Shears

MS-002-28
Sheep Shears
6.5"



Sheep Shears

MS-002-29
Sheep Shears
5"



Sheep Shears

MS-002-30
Sheep Shears -
6 1/2"



Single Bow Sheep Shears

MS-002-31
Sheep Shears
3.5"



Single Bow Sheep Shears

MS-002-32
Sheep Shears

KICK STOPS & COW HOBBLER



**Cow Hobble
MS-002-33**

Cow hobble stops the cow from kicking or switching her tail. Hooks over tendons of hind legs and adjusts to any size cow. Galvanized steel. Wt. 1.5 lb.



**Kow Kant Kick
MS-002-34**

Kow Kant Kick fits any size animal and is now improved with a black powder-coated finish. To use, slip over the animal back and crank the handle to immobilize the muscles of the rear legs that control kicking. It works on either side of the pin bone. The arms are rubber-tipped to protect the animal. Wt. 13 lbs.



**Kick-Stop
MS-002-35**

Kick-Stop has an easy, pushbutton adjustment that allows it to fit any size animal. Slips into place with just a push of the button, and it fits snugly over back for optimum control. Handy and easy to carry, Kick-Stop weighs only 4 lbs. as compared to other anti-kick devices weighing 12 lbs. or more. Scientifically built for long, hard usage, it will be welcomed by every stock owner who handles or treats cattle. Individually boxed with informative label. Wt. 4 lbs.

POLANSKY VAGINAL SPECULUM



Killian Veterinary Vaginal Speculum

MS-002-36

Killian model. Length of jaws 6 cm and 8.5 cm, stainless. 6 cm, 8.5 cm

Veterinary Vaginal Speculum


MS-002-37

Veterinary vaginal speculum 27 cm

Polansky Veterinary Vaginal Speculum

MS-002-38

The speculum has 3 blades. Has butterfly screws with 2 blades on top for wide exposure. 11"



Vaginal Speculum F/Sheep

MS-002-39

Vaginal speculum 20.3 cm

Sheep Vaginal Speculum

MS-002-40

Polansky Vaginal Speculum 7"

HOG CATCHERS & SLAPPERS



**Ring-O-Matic Hog Catcher
MS-002-41**

Made for one-person use, the Ring-O-Matic Hog Catcher has a smooth, sanitary cable on the steel handle. Self-locking feature allows for quick, easy, sure snaring of the animal. Strong enough for large hogs. Wt. 3.5 lbs.



**Hog Holder
MS-002-42**


Hog Holder features 3/16" galvanized aircraft cable, durable high-strength tubing made of nickel-plated rust-resistant carbon steel, and comfortable carbon steel handle. This large Hog Holder is the choice for large hogs. It has a cable stop to avoid kinking, rubber grip for good hand control and a safety lock with quick easy release. Hog Holders are useful for blood testing, giving I.V.s, ringing, tagging or handling in farrowing crates. Wt. 2 lbs. 37"



**Hog Holder
MS-002-43**

Hog Holder features 3/16" galvanized aircraft cable, durable high-strength tubing made of nickel-plated rust-resistant carbon steel, This 16" Hog Holder is the choice for handling hogs in tight quarters. Hog Holders are useful for blood testing, giving I.V.s, ringing, tagging and handling in farrowing rates. Wt. 1 lb. 16"

RESTRAINING EQUIPMENTS



Strub

MS-002-44

Restraining Equipment
30 cm, 1 1/4



Nose Ring

MS-002-45

2" inside diameter, with safety screw
30 cm, 1 1/4



Mouth Gag

MS-002-46

Restraining Equipment
10 1/2 cm, 4 1/8, 15 cm, 6



Restraining Equipment

MS-002-47


Restraining Equipment
23/4



Restraining Equipment

MS-002-48

Restraining Equipment
18 1/2 cm, 7 1/4



Forged Steel

MS-002-49

Restraining Equipment

SELF RETAINING RETRACTORS



Weitlaner
MS-002-50
Sharp/Blunt
13 cm , 5" , 161/2 cm, 61/2



Collin
MS-002-51



Weitlaner
MS-002-53



Weitlaner
MS-002-53
Blunt
13 cm , 5" , 161/2 cm, 61/2



Weitlaner
MS-002-54
Blunt
13 cm , 5" , 161/2 cm, 61 /2



Weitlaner
MS-002-55
Blunt
13 cm , 5" , 161/2 cm, 61 /2

ANIMAL SCISSORS FOR EQUIPS



Cattle-Clipper
MS-002-56
Narrow toothed
17 cm , 6 3/4



Cattle- And Dog Clipper
MS-002-57
Strong model, wide toothed
17 cm , 6 3/4



Doyen Clamp
MS-002-58
Stainless steel
24 cm , 9 1/2



Sheep Shear
MS-002-59
Stainless steel
Large Size 14





Hoof Shears
MS-002-60
Curved on flat scissors
For Sheep, 8



Torsion Fork
MS-002-61
With Fixation bar stainless steel
For Sheep, 8

INSTRUMENTS FOR DOG CARE



Instruments For Dog Care

MS-002-62

Ear clamps for dog, straight of light metal
11 1/2 cm, 4 1/2, 15 1/2 cm, 6 11/8

**Dog Nail Trimmer 17cm blue
handle**

MS-002-63

Stainless Steel

MS-002-64

Instruments for dog Care
17 cm , 6 3/4



Stainless Steel

MS-002-65

Instruments for Dog Care
16 cm, 6 1/4



Stainless Steel

MS-002-66

Curved Scissors 17
cm , 6 3/4



Stainless Steel

MS-002-67

Curved scissors 17
cm , 6 3/4

TOOTH INSTRUMENTS



**Tooth Instrument
MS-002-68**

Tooth forceps for whelps and other small animals, slightly bent stainless steel handle
20 cm, 8



**Flohr
MS-002-69**

Root elevator, right stainless steel handle
20 cm, 8



**Gunther
MS-002-70**

Tooth forceps for large animals
20 cm, 8



**Tooth Rasp
MS-002-71**

With interchangeable blade single piece
55 cm, 22 cm



**Tooth Rasp
MS-002-72**

For lower jaws with interchangeable blade, rasp surface on one side and the file on the other side, Folding (two pieces) straight or bended
55 cm, 22



**Scaler, Hook Shape
MS-002-73**

Tooth Instrument
55 cm, 22


TRAYS FOR INST & SYRINGES



Steriler Drums
MS-002-74



High Polished
MS-002-75



Trays For Inst
MS-002-78
Traysfor Inst
220 X 70 mm



Perforated Trays
MS-002-76
Trays for Inst
18/8, 505 X 305 X 60 mm



Trays
MS-002-77



TRAYS FOR INST & SYRINGES



Glycerine Syringe
MS-002-79 With
3
rings, i cannula,
220 X 70 mm



Pomerory Ears Syringe
MS-002-80
Chrome plated
220 X 70 mm

MANZIM
SURGICO

VARIOUS INSTRUMENTS



Hug

MS-002-81


Teat tumor extactot
15 cm, 6



Paimans

MS-002-82

Teat Dilator



Schecker

MS-002-83



Ullner

MS-002-84

Teat curette
17 cm, 6 3/4



Teat Slitter

MS-002-85

French pattern
10 cm, 4



Teat Dilator With Screw

MS-002-86

Various Instrument
14 1/2 cm, 5 3/4

VARIOUS INSTRUMENTS




Hug Teat Lancet

MS-002-87

Various Instrument

15 1/2 cm, 6 1/8



POST



Virchow
MS-002-88
Stainless Handle
24 1/2cm, 9 1/2



Fig. 1
(180)
Kirkland

Virchow
MS-002-89
Wooden Handle
24 cm , 9 1/2



Fig. 2
(160)
Kirkland

Virchow
MS-002-90
Wooden Handle
25 1/2cm, 9 1/2



Fig. 1
Orson

Post Mortem
MS-002-91
Post mortem
Staniless handle

28 cm, 11"



Charriere
MS-002-92
Adjustable, with 3 saw blades, 3, 5 and
7mm wide 2
28 cm, 11"



Fig. 1/2
TG Orson

Satterlee
MS-002-93
Post Mortem
31 cm, 12 1/4

POST



Fig. 1/2
TG Chisel

Stille

MS-002-94

Post Mortem

23 9"



Fig. 1
TG Chisel

Virchow

MS-002-95

Stainless

22 cm, 8 3/4"

MANZIM
SURGICO

MEASURING INSTRUMENTS



Fig. 2
Orban

Boley Gauge

MS-002-96

Screw lock
10 cm , 4



Vernier Caliper

MS-002-97

Measuring Instrument
11 1/2cm 4 1/2



Veterinary Pal Caliper

MS-002-98

Measuring Instrument
10 cm , 4



Lwanson Clipper

MS-002-99

All uses, for metal- for wax
10 cm , 4



MANZIM SURGICO

We Are Manufacturer & Exporter of Veterinary Instruments.

CONTACT US

ADDRESS: Khatana Street Near Islam Masjid. P/O Box
51310 S-I-E, Shatab Garh, Sialkot, Pakistan.

TEL: 0092-52-3253436

FAX: 0092-52-3253436

MOBILE# 0092-300-5159397
0092-307-5311117

EMAIL: info@manzimsurgico.com
sales@manzimsurgico.com

WEB: www.manzimsurgico.com